

# Mission® HemoPro

## Hemoglobin Microcuvette Package Insert

REF C133-3011-201 # CCS-171 English

For measuring the hemoglobin concentration in capillary or venous whole blood.  
For professional in vitro diagnostic use only.

### INTENDED USE

The Mission® HemoPro Hemoglobin Microcuvette is specifically designed to work with the Mission® HemoPro Hemoglobin Meter to measure the hemoglobin concentration in capillary or venous whole blood. Microcuvettes are designed for professional in vitro diagnostic use only.

### SUMMARY

The major role of hemoglobin is to carry oxygen from the lungs to the tissues and return carbon dioxide from the tissue to the lungs. So, the changes of Hb concentration in the blood can cause several diseases such as anemia and even death. One has to be noted that hemoglobin is a key hematological metric in medical diagnostics, therefore, the accurate determination of this hemo-protein is essential in a number of human pathologies<sup>1</sup>. Hemoglobin assessments are used widely to screen individuals for anemia, to draw inferences about the iron status of populations, and to evaluate responses to nutritional interventions<sup>2</sup>.

The Mission® HemoPro Hemoglobin Testing System can quickly and accurately test the hemoglobin concentration and calculate hematocrit level. The measurement range of the Mission® HemoPro Hemoglobin Testing System is 0.0-26.0 g/dL.

### PRINCIPLE AND REFERENCE VALUES

#### Principle

The system consists of a meter together with microcuvettes. The microcuvette serves as both a pipette and a measuring cuvette which is designed for single use only. A blood specimen or control solution of approximately 10 µL is drawn into the cavity by capillary action. The meter measures whole blood absorbance at the Hb/HbO<sub>2</sub> isobestic point. It measures at two wavelengths (505 and 880 nm) to compensate for turbidity. The Mission® HemoPro testing system is calibrated against the international reference method (HICN), recommended by the ICSH for hemoglobin determination<sup>3</sup>. It requires no further calibration.

#### Reference values

The reference hemoglobin values are listed in the following table:

Unit	Hemoglobin		
	g/dL	g/L	mmol/L
Men	13.0-17.0	130-170	8.1-10.5
Women	12.0-15.0	120-150	7.4-9.3
Children	11.0-14.0	110-140	6.8-8.7

The data in the above table are from Dacie and Lewis Practical Hematology<sup>4</sup>. Reference ranges may vary between laboratories. Every laboratory should establish its own reference range as needed.

### REAGENTS

This product is composed of polystyrene plastic and does not contain any active ingredients.

The performance characteristics of these hemoglobin microcuvettes have been determined in both laboratory and clinical tests. This test has been developed to be specific for the measurement of hemoglobin with the exception of the interferences listed. Refer to the Limitations section for detailed information.

### PRECAUTIONS

- Do not use after expiration date.
- Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal.
- Any serious incident that has occurred in relation to the meter shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### STORAGE AND STABILITY

- The microcuvettes should be stored at 2-40 °C (35.6-104 °F).
- The microcuvettes in the canister can be used until the expiration date printed on the package, regardless of whether the package has been opened or not.
- Always keep the canister closed. All unused microcuvette must be kept in the original package.

### SPECIMEN COLLECTION AND PREPARATION

See User's Manual before using for complete instructions on blood specimen collection.

#### For Venous whole blood specimen:

- EDTA or heparin anticoagulants can be used when collecting the venous whole blood specimen.
- Venous whole blood specimen with EDTA or heparin anticoagulants can be stored in a cool, dry area at 2-30°C (35.6-86°F) for 24 hours. Store them away from heat and direct sunlight.
- Do not freeze. Frozen blood specimen is severely lysed and may lead to incorrect results.
- Some venous whole blood specimens are easily lysed. A red color present in the plasma is a clear indication of blood lysis. If the red color is present, do not use the specimen for testing. Collect a new specimen and conduct the test immediately or within a very short period of time.
- Gently mix the venous whole blood for at least 15 minutes before testing. To avoid blood lysis, do not stir the specimen excessively or mix the specimen for too long.

#### For Capillary blood specimen:

- Prior to testing, wipe the test site with an alcohol swab or wash with soapy water. Then dry hands and test site thoroughly. Any cream or lotion residue on the testing site may affect the test results.
- Do not use iodine tincture disinfectants. The residual iodine tincture on the skin may lead to inaccurate results.
- Capillary blood specimens should be used immediately to avoid blood clotting.

### MATERIALS

#### Materials Provided

- Microcuvettes


- Package Insert

#### Materials Required But Not Provided

- Meter
- Sterile Lancets
- Lancing Device
- Control Solution

### DIRECTIONS FOR USE

Please make sure the meter, microcuvette, venous whole blood specimen, and/or control solution are placed at 10-40°C (50-104 °F) for at least one hour before testing. Refer to User's Manual for detailed instructions.



- A. Filling end
- B. Optical eye
- C. Handheld end

