


NEW

# glucoCcheck

## METER OVERVIEW



MY LIFE, MY DECISIONS.  
**EVOLVE.**

Blood glucose monitoring system



Display Screen

Test Strip Slot  
Insert test strip here to switch the meter on for testing.

Strip Indicator Light  
Light indicates the position of the test strip slot.

Test Strip Ejector  
Eject the used strip by sliding the button upwards.

Main Button (M)  
Enter the meter memory and silence a reminder alarm.

Data Port  
Download test results with a cable connection.

Battery Compartment  
Uses 2x AAA batteries

SET Button (S)  
Enter and confirm the meter settings.

### HOMEMED

Authorised distributor:  
Homemed (Pty) Ltd  
Co Reg 2003/022932/07  
T: 0861 106 150  
Address: 100 Sovereign Drive,  
Route 21 Corporate Park,  
Nellmapius Drive, Irene, 0157,  
South Africa  
W: [www.homemed.co.za](http://www.homemed.co.za)  
E: [sales@homemed.co.za](mailto:sales@homemed.co.za)

EC REP MedNet GmbH  
Borkstraße 10,  
48163 Münster, Germany

TaiDoc Technology Corp.  
B1-7F, No. 127,  
Wugong 2nd Rd., Wugu Dist.,  
24888 New Taipei City, Taiwan



For self-testing  
Model No.: TD-4125



NO CODING

GEQRG ver 1.0 20170420  
315-4125200-004

QUICK REFERENCE GUIDE.  
**EVOLVE.**

# QUICK REFERENCE GUIDE. EVOLVE.

NEW

# glucoCheck

